

Wednesday Standards

5 uL from the 1 ml of incubation medium spotted on a filter

First sample from cells used for pH experiment.

Second sample from cells used for either Val or Glu expts.

Last sample is a blank filter to give the background count level.

145	**	-	1	1	3249.00	3.51	1.00	192.26	62.0	↓ Val
				2	6306.00	2.52				
				3	6311.00	2.52				
146	**	-	2	1	4038.00	3.15	1.00	193.84	61.0	pH/glu
				2	5841.00	2.62				
				3	5849.00	2.62				
147	**	-	3	1	24.00	40.82	1.00	195.42	60.0	BIC
				2	31.00	35.92				
				3	39.00	32.03				